P R A X I S
P R A C T I C E
P R A X I S
P R A C T I C E

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Three-stage fermenter helps Strathmann produce promising vaccines

Strathmann has an eye on the future of the pharmaceutical industry. In its new 5000-square-meter biotechnology center at Dengelsberg, Germany, the company manufactures drugs already on the market as well as biotechnological actives for clinical studies. After the smooth startup of a three-stage fermenter made by Bioengineering, Strathmann is also using it to produce bacterial vaccines against diseases such as AIDS.


The supreme law of pharmaceuticals is that the manufacturing plant, with its highly sensitive process and equipment, must meet the utmost standards of quality and comply with all international regulations such as those of the FDA and EMEA. So Strathmann had good reasons for choosing globally respected specialists with the experience needed to supply and outfit the Dengelsberg center. Erected at a cost of 25 million euro,

the complex brings the former Strathmann Research and Strathmann Biotec laboratories under a single roof.

Plant design for 100-gram scale
The Strathmann strategy: from idea to approved drug. The multi-purpose facility, which houses Bioengineering's three-stage fermenter, makes it possible to produce biotechnological actives and cell-bank systems in various levels of clean rooms. Processes that the customer runs only under lab conditions can be scaled up here. The plant is therefore designed to work on 100-gram scale, and the vaccine production program will soon include work on an AIDS vaccine. The intention is to «outsmart» the virus using a genetic sequence on a cassette holding a mixture of up to 5000 variants of a HIV surface protein.

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Problem-free multiphase start-up
Careful initial startup is a key requirement if a plant is to meet all hygienic and quality standards from the very beginning. After a factory acceptance test (FAT) in the fabrication facilities of Bioengineering AG, the three-stage fermenter began its journey to Germany, where Bioengineering AG carried out a two-week program of mechanical and electrical reassembly at the Dengelsberg complex. The system was brought on-stream in three phases during the fourth quarter of 2002. Commissioning was governed by GMP guidelines and involved Strathmann staff members. The future production team was acquainted with the operating software during the installation qualification (IQ) and functional qualification (FQ) phases.

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Over 120 automated operations
More than 120 fermenter operations are automated and can be initiated with a simple mouse click. During initial startup, all sequences were put through another round of in-depth testing; some minor errors were identified and corrected on the spot. The process centered on the functional design specification (FDS), which had been tailored to individual customer requirements in the design phase. Every control circuit had to be compared with its counterpart in the FDS. Finally, a sterile test run was conducted by Strathmann with Bioengineering staff in attendance. A microbiological assessment showed the fermenter to work correctly at all points, and it was officially handed over to begin the production of promising pharmaceuticals.

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
Fermenter stages: 20, 200 and 1000 liters
 The fermenter built for Strathmann comprises three stages and is designed to grow batch and fed-batch bacterial cultures. The first line is a 20 L fermenter that produces the preculture for the two larger units. The second includes a 200 L fermenter and five feed vessels from which correction media (such as acids, bases, antifoaming agent, ammonia and nutrient medium) can be conveyed. The third unit, a 1000 L fermenter, also has five feed vessels for correction media.

Automated cleaning and sterilization
 The three fermenters are linked by fixed piping. Automated in situ cleaning and sterilization are possible thanks to the permanent connections between tanks and the connections of the feed vessels to the main vessels. In all, the fermenter can execute more than 120 automated operations under computer control.

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